Easycrush 2026 – Product Manual

(Compliant with EU regulations for non-medical devices, HS code 8438.80.9900)

ENGLISH

1. Manufacturer Information

Manufacturer: EasyCrush ApS 📍 Nydamsvej 43, 8362 Hørning, Denmark 📧 info@easycrush.eu 📞 +45 70 20 58 26 / +31 06 38 07 58 84 🌐 www.easycrush.eu

This device is manufactured in compliance with applicable EU directives and regulations.

2. Intended Use

The Easycrush 2026 is an electric pill crusher designed for environments where medication tablets must be crushed safely and efficiently (e.g., care homes, hospitals, private households). ⚠️ Important: • The Easycrush 2026 is not a medical device. • Classified as machinery for preparation of food or pharmaceuticals under HS code 8438.80.9900. • Intended exclusively for crushing solid oral dosage forms (pills, tablets).

3. Safety Instructions

General Safety ⚠️ Use only for intended purpose. 🚫 Do not allow children or unauthorized persons to operate. 🧴 Keep hands, hair, and loose clothing away from moving parts. 👁️ Always supervise during operation.

Electrical Safety ⚡ Use only supplied charger. 💧 Do not immerse in water. ✂️ Damaged cords must be replaced by manufacturer or authorized personnel.

Health & Hygiene 🧴 Always use disposable cups and foil (rondels) to prevent cross-contamination. ✨ Clean exterior with soft, damp cloth; do not use abrasive or aggressive cleaners. ♻️ Do not reuse disposable consumables.

Warnings ⚠️ Do not crush more than 8 pills at a time. 🔒 Do not operate with lid open. 🔥 Stop immediately if unusual sounds, smells, or overheating occur.

4. Technical Specifications

• Product Name: Easycrush 2026 • HS Code: 8438.80.9900 • Cycle Time: ~8 seconds per crushing cycle • Capacity: Small cup 1–4 pills, Large cup 1–8 pills • Motor Power: 2–3× stronger than comparable devices • Power Supply: Rechargeable battery with charger (Input: 100–240V, Output: DC) • Dimensions: [Insert L170mm × W100mm × H194.5mm] • Weight: [2.4 kg] • Noise Level: < [58 dB] • Country of Origin: Denmark

5. Package Contents

📦 Includes: • Easycrush 2026 main unit with two grinding heads • Stainless steel grinding heads (small & large) • Blue safety lid • Charger/power supply • Disposable cups (small & large) • Disposable foil seals (rondels) • User manual

6. Instructions for Use

Step 1 – Preparation 1 Place on stable surface. 2 Press ON/OFF (indicator light on). 3 Remove blue lid (turn counterclockwise).

Step 2 – Loading 4 Select cup size: Small 1–4 pills, Large 1–8 pills.
5 Place pills in disposable cup. 6 Cover cup with foil seal (rondel).
7 Insert grinding head into cup.

Step 3 – Crushing 8 Place prepared cup + rondel + grinding head on base plate. 9 Reattach lid (turn clockwise). 10 Crushing starts automatically. 11 Cycle ~8 seconds. 12 Completion indicated by sound stop & light.

Step 4 – After Use 13 Remove lid, take out cup with grinding head inside. 14 Hold rondel, remove grinding head carefully. 15 Place grinding head back on base. 16 Slide rondel to release powder into cup. 17 Dispose of rondel. 18 Mix powder with water, gel, or food. 19 Administer to patient.

7. Cleaning & Maintenance

20 Clean exterior with soft cloth; use alcohol on grinding head after each use. 21 Do not immerse in water. 22 Store grinding head on base when not in use. 23 Replace disposable cups and rondels after every use. 24 Check lid and grinding head for wear regularly.

8. Storage & Transport

25 Store in dry, cool environment. 26 Avoid sunlight and extreme temperatures. 27 Transport in original packaging.

9. Disposal & Recycling

28 Dispose according to WEEE Directive (2012/19/EU). 29 Do not discard with household waste. 30 Packaging materials are recyclable.

10. Compliance & Conformity

31 Machinery Directive 2006/42/EC 32 Low Voltage Directive 2014/35/EU 33 EMC Directive 2014/30/EU 34 RoHS Directive 2011/65/EU 35 WEEE Directive 2012/19/EU 36 Declaration of Conformity available on request.

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