



DATA SHEET

EASYCRUSH MEDICINE CUP LIDS LARGE & SMALL



**WEBSITE
LINK**

2026

PRODUCT HIGHLIGHTS

EasyCrush disposable medication cup lids are designed for use with the EasyCrush electric pill crusher cups, allowing for optimal storage safety and transfer of medicine.

Pills are placed in the cup, where they are ground into a fine powder by the grinder, ensuring no mixing or contamination occurs between different medications. Once the crushing process is complete, the powder can be combined with a swallowing gel, making it easier to administer with a more enjoyable taste.

TECHNICAL SPECIFICATIONS

LARGE - MADICINE CUP LID ITEM NUMBER: EC-MCLL2000



Material	Low Density Polyethylene, LDPE
Unit Net Weight	1.69 grams
Bundle Net Weight	3.38 kg.
Bundle Gross Weight (incl. packaging)	kg.
Quantity	2000 pcs. in a box
Dimensions	Diameter Ø60mm x Height 4mm
To be used with	EasyCrush Medicine Cup - Large (Item number: EC-MCL2000)
Type	Disposable (one time use)
Certificate	CE
GTIN Code	5745001384067
Country of Origin	China

SMALL - MADICINE CUP LID

ITEM NUMBER: EC-MCLS2000



Material	Low Density Polyethylene, LDPE
Unit Net Weight	1.18 grams
Bundle Net Weight	2.36 kg.
Bundle Gross Weight (incl. packaging)	kg.
Quantity	2000 pcs. in a box
Dimensions	Diameter Ø49mm x Height 4mm
To be used with	EasyCrush Medicine Cup - Small (Item number: EC-MCS2000)
Type	Disposable (one time use)
Certificate	CE
GTIN Code	5745001384081
Country of Origin	China

NOTES

EASYCRUSH MEDICINE CUP LIDS

DESIGNED FOR CARE, BUILT FOR SAFETY

Contact us to learn more!

WEBSITE

www.easycrush.eu

ADDRESS

Birkemosevej 3A,
8361 Hasselager
Denmark

PHONE

T. : +45 123 4567

EMAIL

info@easycrush.eu
